

CRM197 protein, His tag

Cat. No. GTX136284-pro

Applications WB

Package
1 mg

Applications

Observed MW (kDa) 60 kDa.

Properties

Form Lyophilized powder

Buffer Reconstitute with distilled water to 1 mg/ml. Lyophilized from PBS, 5% sucrose.

Preservative No preservatives

Storage Store at 4°C or below. After reconstitution, keep as concentrated solution. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. For long-term storage after reconstitution, aliquot and store at -70°C or below. Do not vortex.

Region/Sequence 26-560aa of CRM197 protein (WP_003850266.1) with N-terminal His tag

Purity >95%

Endotoxin <100 (EU/mg)


Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX136284-pro Image

GTX136284-pro CRM197 Carrier Protein, His tag was analyzed using SDS-PAGE and stained with Coomassie blue and captured by monochrome camera.



For full product information, images and publications, please visit our [website](https://www.genetex.com).



GTX136284-pro WB Image

CRM197 Carrier Protein, His tag (1 μg, GTX136284-pro) was separated by 12% SDS-PAGE, and the membrane was blotted with CRM197 antibody [HL1808] (GTX637460) diluted at 1:3000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).