

SMAD1+5+8(9) antibody

Cat. No. GTX13724

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000

Not tested in other applications.

Calculated MW 52 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.2% Gelatin
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Synthetic peptides corresponding to aa19-33 and 315-329 of human Smad1. Sequence homology: Smad1, Smad5, Smad8, Smad9.
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX13724 WB Image

WB analysis of Saos-2 cells using GTX13724 SMAD1+5+8(9) antibody.



GTX13724 WB Image

WB analysis of Saos-2 cells using GTX13724 SMAD1+5+8(9) antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).