

Bacillus anthracis protective antigen antibody

Cat. No. GTX13809

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Bacillus anthracis

Package
200 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	2 µg/ml - 10 µg/ml

Not tested in other applications.


Properties

Form	Liquid
Buffer	Serum diluted with PBS, 0.05% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	The antibody was developed by immunizing rabbits with synthetic peptides corresponding to amino acids 186-192 of the Anthrax PA protein (approx. 83 kDa).
Purification	Unpurified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX13809 WB Image

WB analysis of recombinant anthrax PA protein using GTX13809 Bacillus anthracis protective antigen antibody.

Lane 1 : pre-bleed

Lane 2 : primary antibody

Dilution : 1:1000



GTX13809 WB Image

WB analysis of recombinant anthrax PA protein using GTX13809 Bacillus anthracis protective antigen antibody.

Lane 1 : pre-bleed

Lane 2 : primary antibody



For full product information, images and publications, please visit our [website](https://www.genetex.com).