

Human IL1 beta protein, His tag

Cat. No. GTX138452-pro
Species Human

Package
100 µg

Applications
Observed MW (kDa) 27 kDa.

Properties
Form Lyophilized powder

Buffer Reconstitute with distilled water to 0.9 mg/ml. Lyophilized from PBS, 2% Trehalose, 5% Mannitol, 0.01% Tween-80.

Preservative No preservatives

Storage Store at 4°C or below. After reconstitution, keep as concentrated solution. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. For long-term storage after reconstitution, aliquot and store at -70°C or below. Do not vortex.

Region/Sequence C-terminal His tagged full-length Human IL1B protein (1-269aa of NP_000567.1)

Expression System HEK293 cells


Purity >95%

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX138452-pro Image

GTX138452-pro Human IL1 beta, His tag was analyzed using SDS-PAGE and stained with coomassie blue and captured by monochrome camera.



For full product information, images and publications, please visit our [website](https://www.genetex.com).