

Human Ferritin Light Chain protein, His tag

Cat. No. GTX138607-pro**Species** Human**Package**
100 µg

Applications

Observed MW (kDa) 26 kDa.


Properties

Form Lyophilized powder**Buffer** Reconstitute with distilled water to 0.3 mg/ml. Lyophilized from PBS, 2% Trehalose, 5% Mannitol, 0.01% Tween-80.**Preservative** No preservatives**Storage** Store at 4°C or below. After reconstitution, keep as concentrated solution. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. For long-term storage after reconstitution, aliquot and store at -70°C or below. Do not vortex.**Region/Sequence** C-terminal His tagged full-length Human Ferritin Light Chain protein (1-175 a.a. of NP_000137.2)**Expression System** HEK293 cells**Purity** >85%**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

**GTX138607-pro Image**

GTX138607-pro Human Ferritin Light Chain protein, His tag was analyzed using SDS-PAGE and stained with coomassie blue and captured by monochrome camera.

For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 10 Page 1 of 1