

Human Prolactin protein, His tag

Cat. No. GTX138611-pro

Applications	ELISA
Species	Human

Package
100 µg

Applications

Observed MW (kDa) 28 kDa.

Properties


Form	Lyophilized powder
Buffer	Reconstitute with distilled water to 1.1 mg/ml. Lyophilized from PBS, 2% Trehalose, 5% Mannitol, 0.01% Tween-80.
Preservative	No preservatives
Storage	Store at 4°C or below. After reconstitution, keep as concentrated solution. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. For long-term storage after reconstitution, aliquot and store at -70°C or below. Do not vortex.
Region/Sequence	C-terminal His tagged Human Prolactin protein (1-227 a.a. of NP_000939.1)
Expression System	HEK293 cells
Purity	>85%
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX138611-pro Image

GTX138611-pro Human Prolactin protein, His tag was analyzed using SDS-PAGE and stained with coomassie blue and captured by monochrome camera.



For full product information, images and publications, please visit our [website](#).



GTX138611-pro ELISA Image

Indirect ELISA analysis was performed by coating the plate with recombinant HEK293 expressed, full-length Human Prolactin protein, His tag (GTX138611-pro) (2.4-153.85 nM). Coated protein was probed with Prolactin antibody (1 µg/mL). Goat anti-rabbit IgG antibody (HRP) (GTX213110-01) (1:10000) and Goat anti-mouse IgG antibody (HRP) (GTX213111-01) (1:10000) were used to detect the bound primary antibody.

Antibodies:

1. Prolactin antibody [GT69] (GTX641658)
2. Prolactin antibody [HL5726] (GTX646677)



For full product information, images and publications, please visit our [website](https://www.genetex.com).