

## IL6 Receptor antibody [15A7]

**Cat. No. GTX14048**

<b>Host</b>	Rat
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2b
<b>Applications</b>	Neutralizing /Inhibition
<b>Reactivity</b>	Mouse

Package  
500 µg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Neutralizing /Inhibition	Assay dependent

**Note : The 15A7 monoclonal antibody has been shown to block the binding of IL-6 to the IL-6 receptor.**

Not tested in other applications.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	Filter-sterilized PBS
<b>Preservative</b>	No preservatives
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	OKT-4 hybridoma cells
<b>Purification</b>	Protein A purified From tissue culture supernatant
<b>Purity</b>	>95% (Determined by SDS-PAGE)
<b>Endotoxin</b>	< 0.002 EU/µg (Determined by LAL assay)
<b>Conjugation</b>	Unconjugated

**Note**  
For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).