

# CD45R antibody [RA3.3A1/6.1 (TIB-146)]

**Cat. No. GTX14090**

<b>Host</b>	Rat
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgM
<b>Application</b>	Depletion
<b>Reactivity</b>	Mouse

**Package**  
500 µg

## APPLICATION

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Depletion	Assay dependent

**Note : The RA3.3A1/6.1 antibody is commonly used for in vivo B cell depletion.**

Not tested in other applications.

## PROPERTIES

<b>Form</b>	Liquid
<b>Buffer</b>	Filter-sterilized PBS
<b>Preservative</b>	No preservative
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Mouse RAW 112 lymphosarcoma cells
<b>Purification</b>	Protein A purified From tissue culture supernatant
<b>Purity</b>	>95% (Determined by SDS-PAGE)
<b>Endotoxin</b>	< 0.002 EU/µg (Determined by LAL assay)
<b>Conjugation</b>	Unconjugated
<b>Note</b>	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.  Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).