

## CD45R antibody [RA3.3A1/6.1 (TIB-146)]

## Cat. No. GTX14090

Host	Rat
Clonality	Monoclonal
Isotype	IgM
Applications	Depletion
Reactivity	Mouse

Package  $500 \, \mu g$ 

## **Applications**

## **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Depletion	Assay dependent
Note: The RA3.3A1/6.1 antibody is comm	nonly used for in vivo B cell depletion.

Not tested in other applications.

Properties	
Form	Liquid
Buffer	Filter-sterilized PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Mouse RAW 112 lymphosarcoma cells
Purification	Protein A purified From tissue culture supernatant
Purity	>95% (Determined by SDS-PAGE)
Endotoxin	< 0.002 EU/
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our website.

Date 2025 / 10 / 28 Page 1 of 1