

Cationic Bovine Serum Albumin (cBSA)

Cat. No. GTX141334-pro

Species Bovine

Package
10 mg

Properties

Form Lyophilized powder

Buffer Reconstitute with 1 ml in distilled water to 10 mg/ml. Lyophilized from distilled water.

Preservative No preservatives

Storage Store at 4°C or below. After reconstitution, keep as concentrated solution. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. For long-term storage after reconstitution, aliquot and store at -70°C or below. Do not vortex.

Purity >95%

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).