

p27 Kip1 antibody

Cat. No. GTX14243

Host	Chicken	
Clonality	Polyclonal	
Isotype	IgY	
Applications	WB	
Reactivity	Human, Mouse, Rat	

Package 50 μg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution		Recommended dilution
WB		Assay dependent
Not tested in other applications.		
Calculated MW	22 kDa. (<u>Note</u>)	

Properties		
Form	Liquid	
Buffer	PBS	
Preservative	0.02% Sodium azide	
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.	
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)	
Immunogen	Immunogen Sequence: GI # 4757962, sequence 1-198 Recombinant cyclin-dependent kinase inhibitor 1B	
Purification	Purified by affinity chromatography	
Conjugation	Unconjugated	
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.	
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.	



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 04 Page 1 of 1

€ 886-3-6208988 📻 886-3-6208989 🐷 infoasia@genetex.com