

M-CSF antibody [5A1]

Cat. No. GTX14316

Host	Rat
Clonality	Monoclonal
Isotype	IgG1
Application	Neutralizing/Inhibition
Reactivity	Mouse

Package
500 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Neutralizing/Inhibition	Assay dependent
Not tested in other applications.	

Calculated MW 61 kDa. ([Note](#))

PROPERTIES

Form	Liquid
Buffer	0.2µm-filtered PBS pH7
Preservative	No preservative
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Purified mouse CSF-1
Purification	Protein G purified From tissue culture supernatant
Purity	>95% (Determined by SDS-PAGE)
Endotoxin	< 0.002 EU/µg (Determined by LAL assay)
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).