

Basal Bodies of Cilia antibody

Cat. No. GTX14373

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC-P, IHC-Fr
Reactivity	Human

Package
100 µg

Applications

Application Note

IHC-P: Use at an assay dependent dilution. IHC-Fr: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the end user.

Product Note This antibody recognises the basal bodies of cilia in all ciliated cells.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.6 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Cytoskeletal preparation of the BHK a21 cell line expressing cilia basal bodies (Human)
Purification	Protein A purified
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).