

Granulocyte antibody [0.N.315]

Cat. No. GTX14497

Host	Mouse
Clonality	Monoclonal
Isotype	lgG1
Applications	ICC/IF, IHC-P, IHC-Fr
Reactivity	Human

Package 250 μl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
Not tested in other applications	

Not tested in other applications.

Product Note Recognizes an unidentified antigen in the cytoplasm of mature granulocytes.

Properties	
Form	Liquid
Buffer	PBS, 0.2% BSA
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Nuclei from pokeweed mitogen-stimulated human peripheral blood lymphocytes.
Purification	Protein G purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 04 Page 1 of 1