

PDE8A antibody

Cat. No. GTX14619

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IP, ELISA, IHC
Reactivity	Human, Mouse, Rat, Bovine, Monkey

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:2000
ICC/IF	1:250
IP	1:250
ELISA	1:10000
IHC	1:250

Not tested in other applications.

Calculated MW 93 kDa. ([Note](#))

Product Note This antibody will label only PDE8A variants. This antibody will not cross react with other PDE family members.

Properties

Form	Liquid
Buffer	Tris/Glycine, 0.5% BSA, 30% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic cyclic peptide from mid-region common to all PDE8A variants.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX14619 WB Image

WB analysis of non-reduced and heated sample lysates using GTX14619 PDE8A antibody.

Lane 1 : SDS-solubilized PDE8 cytosol fraction

Lane 2 : SDS-solubilized treated-PDE8 cytosol fraction

Lane 3 : SDS-solubilized PDE8 cell membrane fraction



For full product information, images and publications, please visit our [website](https://www.genetex.com).