

TCR V gamma 1.1 antibody [2.11]

Cat. No. GTX14669**Host** Armenian Hamster**Clonality** Monoclonal**Isotype** IgG**Applications** FCM**Reactivity** Mouse**Package**

250 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution

FCM

Recommended dilution

Assay dependent

Not tested in other applications.

Product Note

The 2.11 monoclonal antibody reacts with an epitope in the Cr4 domain of TCR Vγ1.1 (T cell receptor V gamma 1.1).

Properties

Form Liquid**Buffer** Filter-sterilized PBS**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.**Concentration** Batch dependent (Please refer to the vial label for the specific concentration.)**Immunogen** 3.13.1 T cell hybridoma**Purification** Protein G purified
From tissue culture supernatant**Purity** >95% (Determined by SDS-PAGE)**Endotoxin** < 0.002 EU/µg (Determined by LAL assay)**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).