

IL33 antibody [Nessy-1]

Cat. No. GTX14709

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, IP, ELISA
Reactivity	Human, Mouse, Rat

Package
50 µg

Applications

Application Note

For IP: 1:100, For WB: 1:1000

Calculated MW	31 kDa. (Note)
Product Note	Recognizes human, mouse and rat pro-IL-33 (precursor) and mature IL-33.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human IL-33 (aa 112-270)
Purification	Protein G-affinity purified.
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).