

Human Thyroglobulin protein

Cat. No. GTX14718-pro

Species Human

References (1)

Package

1 mg

Properties

Form Liquid

Buffer PBS

Preservative No preservatives

Storage Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Expression System Human Thyroid Glands

Purity > 98% by SDS PAGE, ultra-pure, IgG low.

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).