

ADCY4 antibody

Cat. No. GTX14779

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500
ICC/IF	1:100-1:200
IHC-P	Assay dependent
ELISA	1:10000

Not tested in other applications.

Properties

Form	Liquid
Buffer	Tris/Glycine, 0.5% BSA, 30% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.526 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on rat Particulate Adenylyl Cyclase 4 protein.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES



GTX14779 IHC-P Image

IHC-P analysis of mouse colon tissue using GTX14779 ADCY4 antibody.

Dilution : 1:50



GTX14779 IHC-P Image

IHC-P analysis of mouse kidney tissue using GTX14779 ADCY4 antibody.

Dilution : 1:50



GTX14779 WB Image

WB analysis of ADCY4 recombinant protein using GTX14779 ADCY4 antibody.

Dilution : 1:500



For full product information, images and publications, please visit our [website](https://www.genetex.com).