

MERTK antibody

Cat. No. GTX14919

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, FCM, IP, ELISA, IHC
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:750-1:2000
ICC/IF	1:200
FCM	Assay dependent
IP	1:200
ELISA	1:10000
IHC	1:200

Not tested in other applications.

Calculated MW 110 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Tris/Glycine, 0.5% BSA, 30% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.75 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide taken within amino acid region 900-994 on MerTK protein.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTX14919 WB Image

WB analysis of human RPE cell lysate using GTX14919 MERTK antibody.

Dilution : 1:500



GTX14919 IP Image

IP analysis of Long Evans (LE) rat testis tissue lysate using GTX14919 MERTK antibody.

Lane 1 : LE testis extracts immunoprecipitated with primary antibody



For full product information, images and publications, please visit our [website](https://www.genetex.com).