

PICK1 antibody

Cat. No. GTX15001

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Mouse, Rat

Package

50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	200mM Tris, 50mM Glycine, 150mM NaCl, 0.1mM EDTA, 5% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	6xhis tagged fusion protein: EKMELLDQKH VQDIVFLQR FVSTMSKYNN DCYAVLQDAD VFPIEVDLAH TTLAYGPNQG SFTDGEEDE EEEDGAAREV SKDACGATGP TDKGGSWCDS, corresponding to C terminal amino acids 317-416 of Mouse PICK1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).