

Dysferlin antibody (Ready-to-use)

Cat. No. GTX15109

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P, IHC
Reactivity	Human

Package
7 ml

Applications

Application Note

IHC-P: Use at an assay dependent dilution. Perform heat mediated antigen retrieval before commencing with IHC staining protocol. Optimal dilutions/concentrations should be determined by the end user.

Properties


Form	Liquid
Buffer	PBS, 1% BSA
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Immunogen	Synthetic peptide (unfortunately, the amino acid sequence is considered to be commercially sensitive) (Human) (C terminal).
Purification	Unpurified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX15109 IHC-P Image

Human Skeletal Muscle stained with anti-Dysferlin antibody (GTX15108)



For full product information, images and publications, please visit our [website](#).



Human skeletal muscle stained with anti-dysferlin

GTX15109 IHC Image



GTX15109 IHC-P Image

GTX15109 staining Dysferlin in human skeletal muscle by Immunohistochemistry (FFPE-sections).



For full product information, images and publications, please visit our [website](https://www.genetex.com).