

c-Abl antibody

Cat. No. GTX15130

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IP, IHC
Reactivity	Human, Mouse

Package
500 µl

Applications

Application Note

IHC-P: Use at a dilution of 1/50 for 10 min at RT. Boil tissue sections in 10mM citrate buffer, pH 6.0 for 10 min followed by cooling at RT for 20 min. IP: Use at an assay dependent dilution. WB: Use at an assay dependent dilution. Predicted molecular weight: 122.9 kDa. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW 123 kDa. ([Note](#))

Properties


Form	Liquid
Buffer	PBS, 1% BSA
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide from the sequence adjacent to protein tyrosine kinase domain of human c-Abl.
Purification	Immunogen affinity purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

DATA IMAGES



GTX15130 IHC-P Image

GTX15130 staining c-Abl in human placenta by Immunohistochemistry (FFPE-sections).



GTX15130 IHC Image

Human placenta stained with anti-c-Abl



For full product information, images and publications, please visit our [website](#).