

## D4 GDI antibody

**Cat. No. GTX15198**

|                     |            |
|---------------------|------------|
| <b>Host</b>         | Rabbit     |
| <b>Clonality</b>    | Polyclonal |
| <b>Isotype</b>      | IgG        |
| <b>Applications</b> | WB, IHC-P  |
| <b>Reactivity</b>   | Human      |

**Package**  
500 µl

### Applications

#### Application Note

IHC-P: Use at a dilution of 1/100 for 10 min at RT. Boil tissue sections in 10mM citrate buffer, pH 6.0 for 10 min followed by cooling at RT for 20 min. WB: Use at a concentration of 1 µg/ml for 2hours. Predicted molecular weight: 23 kDa. Optimal dilutions/concentrations should be determined by the end user.

**Calculated MW** 23 kDa. ( [Note](#) )

### Properties

|                      |                                                                                                        |
|----------------------|--------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                 |
| <b>Buffer</b>        | PBS, 1% BSA                                                                                            |
| <b>Preservative</b>  | 0.09% Sodium azide                                                                                     |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. |
| <b>Concentration</b> | 0.2 mg/ml (Please refer to the vial label for the specific concentration.)                             |
| <b>Immunogen</b>     | Recombinant full length protein (Human).                                                               |
| <b>Purification</b>  | Immunogen affinity purified                                                                            |
| <b>Conjugation</b>   | Unconjugated                                                                                           |

**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX15198 IHC-P Image

GTX15198 staining D4 GDI in human tonsil by Immunohistochemistry (FFPE-sections).



### GTX15198 IHC-P Image

Human Placenta stained with anti-D4-GDI antibody (GTX15198)



For full product information, images and publications, please visit our [website](#).