

ACE2 peptide (788-805)

Cat. No. GTX15352

Applications	Neutralizing/Inhibition
Species	Human

References (1)

Package

50 µg

Applications

Application Note

For BL: Use at an assay dependent dilution. The peptide is used for blocking the activity of ACE2 antibody. It usually blocks the antibody activity completely in Western blot by incubating the peptide with equal volume of antibody for 30 min at 37°C. Optimal dilutions/concentrations should be determined by the end user.

Properties

Form	Liquid
Buffer	PBS, 0.1% BSA
Preservative	0.02% Sodium azide
Storage	Store the peptide at -20°C, stable for one year.
Concentration	0.20 mg/ml (Please refer to the vial label for the specific concentration.)
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).