

# Human Thyroid (Normal) tissue lysate

**Cat. No. GTX15376**

<b>Species</b>	Human
----------------	-------

<b>Package</b>	100 µg
----------------	--------

**PRODUCT**
**Summary**

Thyroid lysate was prepared by homogenization in lysis buffer. Tissue and cell debris was removed by centrifugation. The product was boiled for 5 min in 1 x SDS sample buffer containing 5% β-mercaptoethanol.

**Tissue type**

Normal

**Properties**
**Form**

Liquid

**Buffer**

1x SDS sample buffer (50mM Tris-HCl, 12.5% Glycerol, 1% Sodium Dodecylsulfate, 0.01% Bromophenol blue), 5% β-mercaptoethanol

**Preservative**

No preservatives

**Storage**

Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

**Concentration**

2 mg/ml (Please refer to the vial label for the specific concentration.)

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).