

PDGF Receptor alpha antibody (Ready-to-use)

Cat. No. GTX15501

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human, Mouse

Package
7 ml

Applications

Application Note

IHC-P: Use at an assay dependent dilution. Staining of formalin fixed tissues requires boiling tissue sections in 10mM citrate buffer, pH 6.0 for 10 min followed by cooling at room temperature for 20 min. Optimal dilutions/concentrations should be determined by the end user.

Properties

Form	Liquid
Buffer	PBS, 1% BSA
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Immunogen	Synthetic peptide (unfortunately, the amino acid sequence is considered to be commercially sensitive) corresponding to C terminal of human PDGFR alpha.
Purification	Unpurified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX15501 IHC-P Image

Human Ovarian Carcinoma stained with anti-PDGFR alpha antibody (GTX15501)


GTX15501 IHC-P Image

GTX15501 staining PDGF Receptor alpha in human ovarian carcinoma by Immunohistochemistry (FFPE-sections).



For full product information, images and publications, please visit our [website](#).