

## Retinoid X Receptor gamma antibody

Cat. No. GTX15518

|              |                   |
|--------------|-------------------|
| Host         | Rabbit            |
| Clonality    | Polyclonal        |
| Isotype      | IgG               |
| Applications | WB, IHC-P         |
| Reactivity   | Human, Mouse, Rat |

References ( 1 )

Package

500 µl

## Applications

## Application Note

IHC-P: Use at a dilution of 1/100. Perform heat mediated antigen retrieval before commencing with IHC staining protocol. WB: Use at a concentration of 1 - 2 µg/ml. Predicted molecular weight: 51 kDa. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW 51 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                             |
|---------------|-----------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                      |
| Buffer        | PBS, 1% BSA                                                                                                                 |
| Preservative  | 0.09% Sodium azide                                                                                                          |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.                      |
| Concentration | 0.2 mg/ml (Please refer to the vial label for the specific concentration.)                                                  |
| Immunogen     | Synthetic peptide (unfortunately, the amino acid sequence is considered to be commercially sensitive) (Mouse) (N terminal). |
| Purification  | Immunogen affinity purified                                                                                                 |
| Conjugation   | Unconjugated                                                                                                                |


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTx15518 IHC-P Image**

Human Prostate Carcinoma stained with anti-Retinoic X Receptor gamma antibody (GTx15518)


**GTx15518 IHC-P Image**

GTx15518 staining Retinoid X Receptor gamma in breast carcinoma by Immunohistochemistry (FFPE-sections).



For full product information, images and publications, please visit our [website](#).