

beta Synuclein antibody

Cat. No. GTX15532

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human, Rat

Package
500 µl

Applications

Application Note

IHC-P: Use at a dilution of 1/100 for 10 min at RT No special pretreatment is required for the immunohistochemical staining of formalin/paraffin tissues Optimal dilutions/concentrations should be determined by the end user.

Product Note This antibody is specific for beta synuclein and does not cross react with alpha or gamma synuclein

Properties


Form	Liquid
Buffer	PBS, 1% BSA
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide (unfortunately, the amino acid sequence is considered to be commercially sensitive) (Human) (C terminal).
Purification	Immunogen affinity purified
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX15532 IHC-P Image

Human Brain stained with anti-Synuclein-beta antibody (GTX15532)



For full product information, images and publications, please visit our [website](#).