

CD11a antibody [121/7]

Cat. No. GTX15818

Host	Rat
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, FCM, IP, Neutralizing /Inhibition
Reactivity	Mouse

Package
250 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
FCM	Assay dependent
IP	Assay dependent
Neutralizing /Inhibition	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.96 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Mouse 180kD alpha chain (CD11alpha) of the integrin Lymphocyte Function-Associated Antigen-1
Purification	Protein G purified
Purity	This product is a Low Endotoxin formulation.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX15818 FCM Image

FACS analysis of mouse splenocytes using GTX15818 CD11a antibody [121/7] compared to an isotype control (blue).

Dilution : 0.5 ug/test incubated for 30 min at 4°C



For full product information, images and publications, please visit our [website](https://www.genetex.com).