

# Cystatin B blocking peptide

**Cat. No. GTX15970-PEP**

<b>Species</b>	Human
----------------	-------

<b>Package</b>	100 µg
----------------	--------

**Applications**

<b>Product Note</b>	The peptide was used in the production of GTX15970.
---------------------	---

**Properties**

<b>Form</b>	Lyophilized powder
-------------	--------------------

<b>Buffer</b>	Reconstitute with 200µl distilled water to obtain a 0.5mg/ml peptide solution. Lyophilized from 5% (v/v) acetonitrile/H <sub>2</sub> O.
---------------	---

<b>Preservative</b>	No preservatives
---------------------	------------------

<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
----------------	--

<b>Region/Sequence</b>	C-QTNKAKHDELT <sub>Y</sub> F
------------------------	------------------------------

<b>Note</b>	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
-------------	---

<b>Note</b>	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
-------------	---



For full product information, images and publications, please visit our [website](#).