

Enterovirus 71 VP1 antibody [1A2-C1]

Cat. No. GTX16114

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Application	ELISA
Reactivity	Enterovirus 71

Package
100 µg

APPLICATION

Product Note The mouse anti-human enterovirus 71-VP1 antibody detects protein derived from HEV without any cross reactivity against unrelated proteins.

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	No preservative
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The immunogen corresponding to HEV 71-VP1 (BAJ05601).
Purification	The IgG fraction of tissue culture supernatant was purified by Protein A/G affinity chromatography.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).