

COX1 antibody [AS70]

Cat. No. GTX16116

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Application	EIA
Reactivity	Human

Package
50 µg

APPLICATION

Application Note

EIA: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW 69 kDa. ([Note](#))

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.08% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Unfortunately, this information is considered to be commercially sensitive
Purification	IgG fraction
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).