

TrpE antibody [17]

Cat. No. GTX16153

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	WB
Reactivity	Species independent

Package
100 µg

Applications

Application Note

WB: Use at a concentration of 1 µg/ml. Detects a band of approximately 47 kDa. Optimal dilutions/concentrations should be determined by the end user.

Product Note GTX16153 does not react with any known mammalian protein.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.08% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Hybridoma produced by the fusion of splenocytes from mice immunized with full length recombinant TrpE protein and mouse myeloma cells.
Purification	IgG fraction
Conjugation	Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX16153 WB Image

TrpE antibody [17]



For full product information, images and publications, please visit our [website](#).