

Galectin 4 antibody [5C8]

Cat. No. GTX16157

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Application	WB
Reactivity	Yeast

Package
50 µg

APPLICATION

Application Note

IP: Use at an assay dependent dilution. WB: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW 99 kDa. ([Note](#))

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.08% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length GAL4 protein expressed in E. Coli.
Purification	Protein G purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).