

Glutathione antibody [D8]

Cat. No. GTX16200

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ICC/IF, IP, ELISA
Reactivity	Species independent

References (1)

Package

100 µg

Applications

Application Note


ELISA, Immunoprecipitation, Western Blot (1:1000, non-reducing), Optimal conditions must be determined individually for each application.

Properties


Form	Liquid
Buffer	PBS
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Purification	Protein A affinity purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX16200 ICC/IF Image

Immunofluorescent analysis of HeLa cells with Glutathione mAb, clone D8, at 1:100 dilution. Red is DAPI nuclear stain.


GTX16200 WB Image

Western blot analysis of Detection of protein-glutathione adducts on Western blot under non-reducing conditions. PTP-1B, phosphotyrosine phosphatase 1B.



For full product information, images and publications, please visit our [website](https://www.genetex.com).