

# DAPI Staining Solution

**Cat. No. GTX16206**

**Applications** WB, ICC/IF, IHC-Fr, FCM, IHC

**References** ( 31 )  
**Package**  
10 ml

## PRODUCT

**Summary** The package of 10 mL is for 100 assays.  
Ex/Em: 340nm / 488nm for DAPI, 360nm / 460nm for DAPI-DNA complexes

## Properties

**Form** Liquid

**Buffer** 8mM Sodium Phosphate, 2mM Potassium phosphate, 140mM NaCl, 10mM KCl.

**Storage** Store at -20°C. Protect from light. When stored as directed, the product should be stable for one year.

**Concentration** 1 µg/ml (Please refer to the vial label for the specific concentration.)

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).