

Goat Anti-Mouse IgG antibody (DyLight488)

Cat. No. GTX16294

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Application	ICC/IF, IHC-P, FACS
Reactivity	Mouse

Package
500 µg

APPLICATION

Application Note

Optimal working dilutions should be determined experimentally by the end user.

PROPERTIES

Form	Liquid
Buffer	PBS, 150mM NaCl, 1% BSA
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Antigen affinity purified
Conjugation	DyLight488
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).