

IL17A antibody [TC11-18H10]

Cat. No. GTX16323

Host	Rat
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, FCM, Neutralizing/Inhibition
Reactivity	Mouse

Package
500 µg

Applications

Calculated MW 17 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	E. coli expressed recombinant mouse IL-17
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).