

Trident 6X Laemmli SDS Sample Buffer

Cat. No. GTX16357

Application WB

Reference (2)
Package
25 ml

PRODUCT

Summary 375 mM Tris-HCl, 9% SDS, 50% glycerol, and 0.03% bromophenol blue.

APPLICATION

Application Note

Add β -mercaptoethanol (9% V/V) before use (reducing).

PROPERTIES

Form Liquid

Storage Store at room temperature.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).