

# Trident 6X Laemmli SDS Sample Buffer

**Cat. No. GTX16357**

**Applications**

WB

**References ( 2 )**

**Package**

25 ml

## PRODUCT

### Summary

375 mM Tris-HCl, 9% SDS, 50% glycerol, and 0.03% bromophenol blue.

### Applications

### Application Note

Add  $\beta$ -mercaptoethanol (9% V/V) before use (reducing).

## Properties

### Form

Liquid

### Storage

Store at room temperature.

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).