

RENT1 antibody

Cat. No. GTX16390

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IP
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
ICC/IF	1:50 - 1:200
IP	1:50 - 1:200

Not tested in other applications.

Calculated MW 124 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide of human UPF1
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 1 of 2

DATA IMAGES




GTX16390 IP Image

IP analysis of HeLa cell lysate using GTX16390 RENT1 antibody.

Antibody amount : 3μg / 300μg lysate

Dilution : 1:1000




GTX16390 ICC/IF Image

ICC/IF analysis of U2OS cells using GTX16390 RENT1 antibody.

Blue : DAPI

Dilution : 1:100



GTX16390 WB Image

WB analysis of various sample lysates using GTX16390 RENT1 antibody.

Dilution : 1:500

Loading : 25μg per lane



For full product information, images and publications, please visit our [website](https://www.genetex.com).