

TTK antibody

Cat. No. GTX16407

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IHC
Reactivity	Human

References (1)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
IHC-P	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Calculated MW	97 kDa. (Note)
Product Note	KO/KD validation is based on published data (PMID: 32121246).

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 517-856 of human TTK (NP_001160163.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX16407 WB Image


WB analysis of HeLa cell lysate using GTX16407 TTK antibody.

Dilution : 1:2000

Loading : 25µg per lane


GTX16407 IHC Image

The data was published in the journal Int J Mol Sci in 2020. [PMID: 32121246](https://pubmed.ncbi.nlm.nih.gov/32121246/)


GTX16407 IHC Image

The data was published in the journal Int J Mol Sci in 2020. [PMID: 32121246](https://pubmed.ncbi.nlm.nih.gov/32121246/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).