

Cystatin C antibody [6F12-C7-D8]

Cat. No. GTX16492

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA
Reactivity	Human

Package
100 µl

Properties


Form	Liquid
Buffer	PBS, 0.5% BSA, 50% glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human cystatin c protein.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX16492 ELISA Image

ELISA analysis of human Cystatin C (4µg/ml) using Cystatin C antibody [6F12-C7-D8] for capture and Cystatin C antibody [7F6-A5-F3] for detection.



For full product information, images and publications, please visit our [website](https://www.genetex.com).