

KCNK1 antibody

Cat. No. GTX16685

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-Fr
Reactivity	Human, Mouse, Rat

Package

50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-Fr	Assay dependent

Not tested in other applications.

Calculated MW 38 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA
Preservative	0.025% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.8 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide RQELRKLKRRFLEEHEC, corresponding to amino acid residues 53-69 (Extracellular, near the P1 loop) of human KCNK1 (Accession : O00180).
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX16685 WB Image


WB analysis of rat brain lysate using GTX16685 KCNK1 antibody preincubated with or without immunogen peptide.

Dilution : 1:200


GTX16685 WB Image

WB analysis of HEK-293-KCNK1 transfected cell lysates using GTX16685 KCNK1 antibody preincubated with or without immunogen peptide.

Dilution : 1:300


GTX16685 IHC-Fr Image

IHC-Fr analysis of mouse cerebellum tissue using GTX16685 KCNK1 antibody.

Panel A : KCNK1 channel appears in glial processes (red).

Panel B : Staining of Purkinje nerve cells with mouse anti-calbindin D28K (a calcium binding protein, green).

Panel C : Merge of KCNK1 channel and calbindin D28K demonstrates the separate localization of these proteins.



For full product information, images and publications, please visit our [website](https://www.genetex.com).