

BRCA1 antibody [MS110]

Cat. No. GTX16780

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-P, IHC-Fr, IP
Reactivity	Human

References (1)

Package

50 µg

Applications

Application Note

For IHC-Fr, IHC-P, IP, IF: Use at an assay dependent dilution. For WB: Use at a concentration of 1.5 - 2.5 µg/ml. Predicted molecular weight: 207 kDa. Optimal dilutions/concentrations should be determined by the researcher.

Calculated MW 208 kDa. ([Note](#))**Product Note** GTX16780 doesnot cross-react with blood group antigens or growth factor receptors such as EGFR.

Properties

Form	Liquid
Buffer	50mM Sodium Phosphate, 0.2% Gelatin
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant full length protein (Human).
Purification	From hybridoma culture supernatant
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).