

## HA tag antibody [12CA5]

Cat. No. GTX16918

|              |                     |
|--------------|---------------------|
| Host         | Mouse               |
| Clonality    | Monoclonal          |
| Isotype      | IgG2b               |
| Applications | WB, IP              |
| Reactivity   | Species independent |

References ( 4 )

Package

200 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:1000               |
| IP                 | Assay dependent      |

Not tested in other applications.

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS                                                                                                                                                                                                                        |
| Preservative  | 0.09% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| Immunogen     | KLH conjugated synthetic peptide encoding HA tag (YPYDVPDYA) was used as antigen.                                                                                                                                          |
| Purification  | Protein G purified, followed by dialysis against PBS.                                                                                                                                                                      |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX16918 WB Image

WB analysis of HA-tag recombinant protein using GTX16918 HA tag antibody [12CA5].



### GTX16918 WB Image

WB analysis of HA-tagged recombinant protein using GTX16918 HA tag antibody [12CA5].



For full product information, images and publications, please visit our [website](#).