

AIP1 antibody

Cat. No. GTX17024

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	2.5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 44 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Aip1 antibody was raised against a 17 amino acid synthetic peptide near the carboxy terminus of the human Aip1. The immunogen is located within amino acids 290 - 340 of Aip1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).


DATA IMAGES

**GTX17024 IHC-P Image**

IHC-P analysis of human brain tissue using GTX17024 AIPL1 antibody.
Working concentration : 20 µg/ml

**GTX17024 IHC-P Image**

IHC-P analysis of human brain tissue using GTX17024 AIPL1 antibody.
Working concentration : 2.5 µg/ml

**GTX17024 WB Image**

WB analysis of human brain tissue lysate using GTX17024 AIPL1 antibody.
Working concentration : (A) 1 and (B) 2 µg/ml



For full product information, images and publications, please visit our [website](#).