

Nkx3.2 antibody

Cat. No. GTX17098

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 35 kDa. ([Note](#))

Product Note This antibody is predicted not to cross-react with other NKX homeobox proteins.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Nkx3.2 antibody was raised against a 19 amino acid synthetic peptide near the center of human Nkx3.2. The immunogen is located within amino acids 180 - 230 of Nkx3.2.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx17098 WB Image


WB analysis of human brain tissue lysate in (A) the absence and (B) the presence of blocking peptide using GTx17098 Nkx3.2 antibody.

Working concentration : 1 µg/ml


GTx17098 IHC-P Image

IHC-P analysis of mouse brain tissue using GTx17098 Nkx3.2 antibody.

Working concentration : 20 µg/ml


GTx17098 IHC-P Image

IHC-P analysis of mouse brain tissue using GTx17098 Nkx3.2 antibody.

Working concentration : 5 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).