

AP3B2 antibody

Cat. No. GTX17145

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 119 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	AP3B2 antibody was raised against an 18 amino acid synthetic peptide near the amino terminus of human AP3B2. The immunogen is located within the first 50 amino acids of AP3B2.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX17145 WB Image

WB analysis of rat brain tissue lysate using GTX17145 AP3B2 antibody.
Working concentration : (A) 1 and (B) 2 μ g/ml




GTX17145 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX17145 AP3B2 antibody.
Working concentration : 20 μ g/ml



GTX17145 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX17145 AP3B2 antibody.
Working concentration : 5 μ g/ml



GTX17145 WB Image

WB analysis of mouse brain tissue lysate using GTX17145 AP3B2 antibody.
Working concentration : (A) 1 and (B) 2 μ g/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).